



Centrifuge Bundle Guide



Simplify your centrifuge selection process with the OHAUS bundle guide! Featuring combinations designed around specific Research or Clinical applications, OHAUS has the right packages to help you **Do More** in your laboratory.

Ingenuously Practical

Research Applications

High-Speed Micro Centrifuge Bundles



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30130866	Centrifuge Micro, 230V, FC5515	1	15 200	19 165×g	24
30642371	Rotor Angle 24×1,5/2,0ml BIOSEALS, V2	1			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30130868	Centrifuge Micro, 230V, FC5515R	1	15 200	19 165×g	24
30642371	Rotor Angle 24×1,5/2,0ml BIOSEALS, V2	1			

Hematocrit Bundle



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30393187	Centrifuge Micro, 230V, FC5513	1	13 500	16 085×g	24 Capillaries
30472304	Rotor, Hematocrit, 24×Capillaries, V2	1			



Microplate Bundle



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314810	Centrifuge Multi, 230V, FC5714	1	4500	2716×g	6 Standard Microplates or 2 Deepwell Plates (max. height 63mm)
30314824	Swing Out Rotor 2×3 MTP with Bucket	1			


Low Volume – Cell Culture Bundles



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314812	Centrifuge Multi, 230V, FC5718	1	5000	3885×g	
30314823	Swing Out Rotor, 4×100 ml Sealable	1			
30314882	Bucket, 3×15 ml D17 mm, conical bottom, w/Rack, 2/pk	2			12
30314881	Bucket, 2×50 ml D29 mm, conical bottom, w/Rack, 2/pk	2			8



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314814	Centrifuge Multi, 230V, FC5718R	1	5000	3885×g	
30314823	Swing Out Rotor, 4×100 ml Sealable	1			
30314882	Bucket, 3×15 ml D17 mm, conical bottom, w/Rack, 2/pk	2			12
30314881	Bucket, 2×50 ml D29 mm, conical bottom, w/Rack, 2/pk	2			8



Mid-Volume – Cell Culture Bundles



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314816	Centrifuge Multi, 230V, FC5816	1	4500	37 355×g	
30314828	Rotor Swing Out, 4×250 ml, ID	1			
30314583	Rack, 10×15 ml D17 mm, conical bottom, 2/ pk	2			40
30304375	Rack, 4×50 ml D29 mm, conical bottom, 2/ pk	2			16



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314818	Centrifuge Multi, 230V, FC5816R	1	4500	37 355×g	
30314828	Swing Out Rotor, 4×250 ml, ID	1			
30314583	Rack, 10×15 ml D17 mm, conical bottom, 2/ pk	2			40
30304375	Rack, 4×50 ml D29 mm, conical bottom, 2/ pk	2			16

High Volume – Cell Culture Bundles



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30553036	Centrifuge Multi, 230V, FC5916	1	4000		
30553085	Swing Out Rotor, 4×750 ml ID Sealable	1			
30602502	Bucket 1 × 750 ml D99 mm, flat bottom, 2/ pk, V1	2			
30553130	Rack 7×50 ml D29 mm, conical bottom, 2/ pk	2		3 416×g	28
30553134	Rack 17×15 ml D17 mm, conical bottom, 2/ pk	2		3 398×g	68



Optional Lid Assembly 30553122 for bucket 30602502 (2/ pk)







Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30553101	Centrifuge Multi, 230V, FC5916R	1	4500		
30553085	Rotor Swing Out, 4 x 750 ml ID Sealable	1			
30602502	Bucket 1 x 750 ml D99 mm, flat bottom, 2/ pk, V1	2			
30553130	Rack 7 x 50 ml D29 mm, conical bottom, 2/ pk	2		3324 x g	28
30553134	Rack 17 x 15 ml D17 mm, conical bottom, 2/ pk	2		4301 x g	68

Optional Lid Assembly 30553122 for bucket 30602502 (2/ pk)





Clinical Applications

Low Volume – Clinical Bundles (Angle Rotor)



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30332131	Centrifuge, 230V, FC5706 with Angle Rotor, 12×15 ml, round bottom/conical bottom	1	6 000	4 265×g	12
30130889	Adapter, D13,5 mm, 7 ml Tubes, 2/pk	6			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30332131	Centrifuge, 230V, FC5706 with Angle Rotor, 12×15 ml, round bottom/conical bottom	1	6 000	3 621×g	12
30130890	Adapter, D13,5 mm, 5 ml Tubes, 2/pk	6			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30393189	Centrifuge Multi, 230V, FC5707 with Angle Rotor 8×15 ml	1	6 800	4 290×g	8
30130889	Adapter, D13,5 mm, 7 ml Tubes, 2/pk	4			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30393189	Centrifuge Multi, 230V, FC5707 with Angle Rotor 8×15 ml	1	6 800	3 618×g	8
30130890	Adapter, D13,5 mm, 5 ml Tube, 2/pk	4			



Low Volume – Clinical Bundles (Swing Out Rotor)



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314810	Centrifuge Multi, 230V, FC5714	1	4500	3214×g	28
30314822	Swing Out Rotor, 4×145 ml, ID	1			
30314857	Rack, 7×5-7 ml D13 mm, Vacutainer®, 2/pk	2			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314810	Centrifuge Multi, 230V, FC5714	1	4500	3237×g	20
30314822	Swing Out Rotor, 4×145 ml, ID	1			
30314858	Rack, 5×4-10 ml D16 mm, Sarstedt®, 2/pk	2			

Mid Volume – Clinical Bundles (Swing Out Rotor)



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314814	Centrifuge Multi, 230V, FC5718R	1	4500	3237×g	20
30314822	Swing Out Rotor, 4×145 ml, ID	1			
30314858	Rack, 5×4-10 ml D16 mm, Sarstedt®, 2/pk	2			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314814	Centrifuge Multi, 230V, FC5718R	1	4500	3237×g	20
30314822	Swing Out Rotor, 4×145 ml, ID	1			
30314857	Rack, 7×5-7 ml D13 mm, Vacutainer®, 2/pk	2			





Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314814	Centrifuge Multi, 230V, FC5718R	1	5 000	4 080×g	32
83041011	Swing Out Rotor, 4×200 ml ID Sealable	1			
83041013	Bucket, 200 ml w/ Cap BIOSEALS, 2/pk	2			
83041028	Rack, 8×7 ml D13 mm, round bottom/Vacutainer®, 2/pk	2			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314814	Centrifuge Multi, 230V, FC5718R	1	5 000	4 080×g	32
83041011	Swing Out Rotor, 4×200 ml ID Sealable	1			
83041013	Bucket, 200 ml w/ Cap BIOSEALS, 2/pk	2			
83041027	Rack, 7×15 ml D17 mm, round bottom/Vacutainer®, 2/pk, V1	2			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314816	Centrifuge Multi, 230V, FC5816	1	4 500	3 689×g	32
30314820	Swing Out Rotor, 4×250 ml ID	1			
30304367	Rack, 8×4-10 ml D16 mm, Vacutainer®, 2/pk	2			



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30314818	Centrifuge Multi, 230V, FC5816R	1	4 500	3 689×g	32
30314820	Swing Out Rotor, 4×250 ml ID	1			
30304367	Rack, 8×4-10 ml D16 mm, Vacutainer®, 2/pk	2			



High Volume – Clinical Bundles (Swing Out Rotor)



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30553036	Centrifuge Multi, 230V, FC5916	1	4000	3 345×g	108
30553085	Swing Out Rotor, 4×750 ml ID Sealable	1			
30602502	Bucket, 1×750 ml D99 mm, flat bottom, 2/pk, V1	2			
30553139	Rack, 27×1,8/7 ml D13 mm, Vacutainer®/Sarstedt®, 2/pk	2			

Optional Lid Assembly 30553122 for bucket 30602502 (2 /pk)




Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30553101	Centrifuge Multi, 230V, FC5916R	1	4500	4 233×g	108
30553085	Swing Out Rotor, 4×750 ml ID Sealable	1			
30602502	Bucket, 1×750 ml D99 mm, flat bottom, 2/pk, V1	2			
30553139	Rack, 27×1,8/7 ml D13 mm, Vacutainer®/Sarstedt®, 2/pk	2			

Optional Lid Assembly 30553122 for bucket 30602502 (2 /pk)



Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30553036	Centrifuge Multi, 230V, FC5916	1	4000	3 398×g	72
30553085	Swing Out Rotor, 4×750 ml ID Sealable	1			
30602502	Bucket, 1×750 ml D99 mm, flat bottom, 2/pk, V1	2			
30553140	Rack, 18×4/10 ml D16 mm, Vacutainer®/Sarstedt®, 2/pk	2			

Optional Lid Assembly 30553122 for bucket 30602502 (2 /pk)





Item N°	Material	QTY	RPM	RCF	Tubes per rotor
30553101	Centrifuge Multi, 230V, FC5916R	1	4500	4301 ×g	72
30553085	Swing Out Rotor, 4×750 ml ID Sealable	1			
30602502	Bucket, 1×750 ml D99mm, flat bottom, 2/pk, V1	2			
30553140	Rack, 18×4/10 ml D16 mm, Vacutainer®/Sarstedt®, 2/pk	2			

Optional Lid Assembly 30553122 for bucket 30602502 (2 / pk)





OHAUS Corporation

Headquartered in Parsippany, New Jersey, OHAUS Corporation is a global leader in manufacturing an extensive line of weighing products, laboratory equipment and analytical instruments that meet and surpass the weighing and measurement needs of a broad range of industries, including laboratory, industrial, educational, food preparation and retail markets. An ISO 9001:2015 manufacturer, OHAUS develops products that are precise, reliable, affordable and backed by industry-leading customer support.

OHAUS Europe GmbH

e-mail: ssc@ohaus.com

Tel: 0041 22 567 53 19

e-mail: tsc@ohaus.com

Tel: 0041 22 567 53 20

www.ohaus.com

*The management system governing
the manufacture of this product
is ISO 9001:2015 certified.*

Ingeniously Practical

80776912_A 20240325
© Copyright OHAUS Corporation